


Product Datasheet

TAB3 Antibody (orb248169)

Description	TAB3, a TAB2-like protein, is a mediator of TAK1 activation in IL-1 and TNF signaling. IL-1 and TNF ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the N-terminal region of human TAB3 was used as the immunogen for this antibody.
Preservatives	1X PBS, pH 7.4
Storage	Aliquot the TAB3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Provided assay concentrations are suggestions only, antibody titration may be required for optimal results.
Formula	0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	Protein A purified antibody
Hazard Information	This TAB3 antibody is available for research use only.
Entrez	257397
Dilution Range	Western blot: 4-8ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human samples u...