

Product Datasheet

Cyclin A2 Antibody / CCNA2 (orb248207)

Description

This antibody recognizes a protein of 50-55

epitope is loc...

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P

Applications
Immunogen

Full length bovine protein was used as the immunogen for this CCNA2 antibody. The epitope has been localized to amino acids 144-148 of the human protein.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Cyclin A2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cyclin A2 antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This Cyclin A2 antibody is available for research use only.

Entrez

890

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.