

Product Datasheet

HER2 Antibody / ErbB2 (orb248329)

Description

This mAb is specific to HER2 / ErbB2 and shows


Species/Host

fami...

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

ELISA, FACS, IF

Applications


Immunohistochemical
analysis of formalin...

Immunogen

Recombinant human protein was used as the immunogen for this HER2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage


Store the HER2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the mAb to be titrated up or down for optimal performance.



Immunohistochemical
analysis of formalin...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This HER2 antibody is available for research use only.

Entrez

[2064](#)

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-4ug/ml

Expiration Date

12 months from date of receipt.