


Product Datasheet

Estrogen Receptor alpha antibody (PE) (orb248340)

Description

Mouse monoclonal antibody against Estrogen

Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FACS
Immunogen	Recombinant human Estrogen Receptor alpha protein (aa 2-185) was used as the immunogen for this antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Estrogen receptor antibody at 2-8°C. Conjugate is light sensitive, store in the dark.
Note	For research use only
Application notes	The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Hazard Information	This Estrogen receptor antibody is available for research use only.
Entrez	2099
Dilution Range	Flow Cytometry: 5ul/test/million cells in 100ul or 5ul/test/100ul of whole blood
Expiration Date	12 months from date of receipt.



FACS
staining of
human
MCF-7
using
Estro...