


Product Datasheet

Estrogen Receptor beta antibody (PE) (orb248345)

Description

Mouse monoclonal antibody against Estrogen

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FACS
Immunogen	The C-terminus fragment of recombinant human Estrogen Receptor beta protein was used as the immunogen for this antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Estrogen Receptor beta antibody at 2-8°C. Conjugate is light sensitive, store in the dark.
Note	For research use only
Application notes	The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Estrogen Receptor beta antibody to be titrated up or down for optimal performance.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	Q92731
Hazard Information	This Estrogen Receptor beta antibody is available for research use only.
Entrez	2100
Dilution Range	Flow Cytometry: 5ul/test/million cells or 5ul/test/100ul of whole blood, Immunofluorescence: 1:50-1:100 for 30 minutes at RT (1)
Expiration Date	12 months from date of receipt.



FACS staining of MCF-7 cells using Estro...