

Product Datasheet

**VEGF Receptor 1 / VEGFR1 Antibody / FLT1
(orb248355)**

Description

VEGFR1 is high affinity receptors for

(VEGFs). On the basis of s...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, Functional Studies, IF

Applications

Immunogen

Recombinant human protein was used as the immunogen for this antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the VEGFR1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Variations in protocols, secondaries and substrates may require the antibody to be titrated for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This VEGFR1 antibody is available for research use only.

Entrez

[2321](#)

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Functional studies (order BSA/sodium azide-free format)

Expiration Date

12 months from date of receipt.