

Product Datasheet

HSP60 Antibody (orb248382)

Description

This antibody recognizes a 60kDa protein, identified as the



Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

IF, IHC-P, WB

Western blot analysis of HeLa cell lysat...

Immunogen

Recombinant human HSP60 protein was used as the immunogen for this antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the HSP60 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This HSP60 antibody is available for research use only.

Entrez

[3329](#)

Dilution Range

Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.