


Product Datasheet

CD11a Antibody / LFA-1 (orb248414)

Description

This antibody recognizes a protein of 180kDa, identified



FACS staining
of human
PBMCs using
CD11a...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IF

Applications
Immunogen

Stimulated human leukocytes were used as the immunogen for this CD11b antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CD11a antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This CD11a antibody is available for research use only.

Entrez

[3683](#)

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml Functional studies (order BSA/sodium azide-free format)

Expiration Date

12 months from date of receipt.