



Product Datasheet

LPAR4 antibody (orb24867)


Description	LPAR4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Lysophosphatidic acid receptor 4 protein (1-43AA)
Target	LPAR4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99677
Dilution Range	WB:1:1000-5000, IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse kidney ti...



Immunohistochemical staining of human en...



Immunofluorescence analysis of HepG2 cel...