


Product Datasheet

Integrin beta 1 antibody (FITC) [Out of stock] (orb248745)

Description	Mouse monoclonal antibody against Integrin beta 1 conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FACS, ICC, IF
Immunogen	Native purified human CD29.
Target	ITGB1
Preservatives	IgG1. Liquid in PBS, pH 7.3, 0.2% rAlbumin, and 0.02% sodium azide.
Storage	Shipped and store at 4°C for one year. Do not freeze.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P05556
Entrez	3688
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
human
periph...