

Product Datasheet

LRRTM2 antibody (orb24887)

Description	LRRTM2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Leucine-rich repeat transmembrane neuronal protein 2 protein (165-422AA)
Target	LRRTM2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43300
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.