


Product Datasheet

Human Ig light chain Kappa antibody [Out of stock] (orb248957)

Description	Mouse monoclonal antibody to Human Ig light chain Kappa
Species/Host	Mouse
Tested Applications	FACS, ICC, IF
Immunogen	Native purified human Ig light chain Kappa.
Preservatives	IgG1. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
human
periph...