

Product Datasheet

LY6G6C antibody (orb24914)

Description	LY6G6C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Lymphocyte antigen 6 complex locus protein G6c protein (19-99AA)
Target	LY6G6C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95867
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.