

Product Datasheet

MAP3K12 antibody (orb24942)

Description	MAP3K12 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 12 protein (101-400AA)
Target	MAP3K12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12852
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.