



Product Datasheet

MAP3K8 antibody (orb24967)


Description	MAP3K8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 8 protein (1-300AA)
Target	MAP3K8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41279
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human glioma tissue.



Immunofluorescence analysis of HepG2 cells.