

Product Datasheet

MPZL2 antibody (FITC) (orb25050)

Description	MPZL2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Myelin protein zero-like protein 2 protein (30-140AA)
Target	MPZL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60487
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.