



Product Datasheet

NDUFB6 antibody (orb25099)

Description	NDUFB6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 protein (2-67AA)
Target	NDUFB6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95139
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of HEK293 W...



Immunohistochemical staining of human lu...