


## Product Datasheet

### AGN 194310 (orb251371)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | AGN 194310(VTP-194310) is a potent and selective pan-RARs agonist with Kd values of 3/2/5 nM for RAR... |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Formula</b>         | C <sub>28</sub> H <sub>24</sub> O <sub>2</sub> S                                                        |
| <b>MW</b>              | 424.55                                                                                                  |
| <b>CAS Number</b>      | [229961-45-9]                                                                                           |
| <b>SMILES</b>          | CC(S1)(C)C=C(C2=CC=C(CC)C=C2)C3=C1C=CC(C#CC4=CC=C(C(O)=O)C=C4)=C3                                       |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                         |



Chemical  
structure of  
AGN 194310