

## Product Datasheet

### IKK(epsilon)-IN-1 (IKKE-IN-1) (orb251384)

|                        |                                                                                                                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | IKKε-IN-1 is a potent IKKε inhibitor; inhibit the in-situ IKK ε-mediated phosphorylation of IRF...                                                                                                           |
| <b>Note</b>            | For research use only                                                                                                                                                                                        |
| <b>Formula</b>         | C <sub>26</sub> H <sub>27</sub> N <sub>5</sub> O <sub>3</sub>                                                                                                                                                |
| <b>MW</b>              | 457.5243                                                                                                                                                                                                     |
| <b>CAS Number</b>      | [1292310-49-6]                                                                                                                                                                                               |
| <b>SMILES</b>          | <chem>O=C([H])([H])C([H])([H])C([H])(C([H])([H])C1([H])[H])OC1C([H])=C([H])C(C2C([H])=C([H])N=C(N=2)N([H])C2C([H])=C([H])C(=C([H])C=2[H])N2C([H])([H])C([H])([H])OC([H])([H])C2([H])[H])=C([H])C=1C#N</chem> |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                              |