

Product Datasheet

HMGCR (phospho-Ser872) antibody (orb251472)

Description

HMGR (phospho-Ser872) antibody

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human HMGR (phospho-Ser872). Please contact us for the exact immunogen sequence. The peptide is available as orb374856.

Target

HMGR (phospho-Ser872)

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

97 kDa

Uniprot ID

[P04035](#), [Q09137](#), [Q5EG47](#)

NCBI

[11](#), [21](#), [000859](#), [000850](#)

Entrez


[3156](#)

Dilution Range


IHC-P: 1:30-100, WB: 1:200

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of paraffin...