



Product Datasheet

BMP2 antibody (orb251474)


Description	BMP2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human BMP2. Please contact us for the exact immunogen sequence. The peptide is available as orb374858.
Target	BMP2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	13 kDa
Uniprot ID	P12643 , P21274 , P49001
NCBI	21 , 001200 , 001191 , 11
Entrez	650
Dilution Range	WB: 1:200, IHC-P: 1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



IF analysis of mouse lung tissue using B...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat sple...