

Product Datasheet

IL11 antibody (orb251476)

Description

IL11 antibody

Species/Host

Rabbit

Reactivity

Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

Applications

ELISA, ICC, IF, IHC-P

Immunogen

KLH conjugated synthetic peptide derived from human IL11. Please contact us for the exact immunogen sequence. The peptide is available as orb374860.

Target

IL11

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

19 kDa

Uniprot ID

[P47873](#), [Q99MF5](#), [P20809](#)

NCBI

[000641](#), [000632](#), [31](#), [11](#)

Entrez


[3589](#)

Dilution Range


IHC-P: 1:200, IF/ICC: 1:200

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of rat subcutaneous tissue t...



IF analysis of rat skin tissue using IL1...