

Product Datasheet

Islet 1 antibody (orb251477)

Description

Islet 1 antibody

Species/Host

Rabbit

Reactivity

Equine, Gallus, Human, Mouse, Rat, Sheep

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Islet 1. Please contact us for the exact immunogen sequence. The peptide is available as orb374861.

Target

Islet 1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

39 kDa

Uniprot ID

[P61372](#), [P61371](#), [P61374](#)

NCBI

[002193](#), [21](#), [002202](#)

Entrez


[3670](#)

Dilution Range


WB: 1:200-1000, IHC-P: 1:200

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat panc...



Western blot analysis of Rat liver tissu...



IHC-P image of rat pancreas tissue using...