

Product Datasheet

CD19 antibody (orb251482)

Description

CD19 antibody

Species/Host

Rabbit

Reactivity

Bovine, Equine, Guinea pig, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, IHC-P

Immunogen

KLH conjugated synthetic peptide derived from human CD19. Please contact us for the exact immunogen sequence. The peptide is available as orb374866.

Target

CD19

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW



61 kDa

Uniprot ID
P15391, P25918, Q5FVF7
NCBI
51, 001761, 001770, 31
Entrez
930
Dilution Range


IHC-P: 1:100, IF/ICC: 1:100

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of rat spleen...



Immunohistochemical staining of paraffin...



IF analysis of rat spleen tissue using a...