

Product Datasheet

Hepatocyte growth factor HGF Antibody (orb251503)

Description

Hepatocyte growth factor HGF Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of mouse HGF (33-64aa QKKRRNTLHEFKKSAKTTLTKEPLLKIKTKK), different from the related human sequence by five amino acids, and from the related rat sequence by one amino acid.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for HGF is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

82945 MW

Uniprot ID


[Q08048](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Mouse

Expiration Date

12 months from date of receipt.



Western blot analysis of NIH Cell using ...