

Product Datasheet

Transcription factor MafA MAFA Antibody (orb251505)

Description

Transcription factor MafA MAFA Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human MAFA (136-167aa EDAVEALIGSGHHGAHHGAHPAAAAAYEAFR), identical to the related mouse sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for MAFA is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

36982 MW

Uniprot ID


[Q8NHW3](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Rat

Expiration Date

12 months from date of receipt.



WB analysis of MAFA using anti-MAFAantib...