

Product Datasheet

Synaptopodin/SYNPO Antibody (orb251508)

Description

Synaptopodin/SYNPO Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Synaptopodin (629-658aa EKPKVTPNPDLLDLVQTADKRRQRDQGEV), different from the related mouse sequence by one amino acid, and identical to the related rat sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for Synaptopodin is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

99463 MW

Uniprot ID


[Q8N3V7](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse

Expiration Date

12 months from date of receipt.



Western blot analysis using Anti-Synapto...