

Product Datasheet

RBP4 Antibody (orb251543)

Description

RBP4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E.coli-derived human RBP4 recombinant protein (Position: E19-201L). Human RBP4 shares 86% amino acid (aa) sequence identity with both mouse and rat RBP4.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for RBP4 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

23010 MW

Uniprot ID


P02753

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat, - Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, By Heat ELISA , 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis using Anti-RBP4 Pi...



IHC analysis of RBP4 using anti-RBP4 ant...



IHC analysis of RBP4 using anti-RBP4 ant...