

Product Datasheet

SFRP2 Antibody (orb251545)

Description	SFRP2 Antibody	130KD — 70KD — 55KD — 35KD — 25KD — 15KD — 1 2
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	E.coli-derived human SFRP2 recombinant protein (Position: L104-C295). Human SFRP2 shares 99% amino acid (aa) sequence identity with mouse SFRP2.	Western blot analysis using Anti-SFRP2 P...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	WB: The detection limit for SFRP2 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	33490 MW	
Uniprot ID	Q96HF1	
Dilution Range	Western blot, 0.1-0.5µg/ml, Human, -ELISA , 0.1-0.5µg/ml, Human	
Expiration Date	12 months from date of receipt.	