


Product Datasheet

LRIG1 antibody (orb251608)

Description	LRIG1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human LRIG1. Please contact us for the exact immunogen sequence. The peptide is available as orb374868.
Target	LRIG1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	120 kDa
Uniprot ID	Q96JA1 , P70193
Entrez	26018
Dilution Range	IHC-P: 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.