

Product Datasheet

Thrombospondin 2 antibody (orb251609)

Description

Thrombospondin 2 antibody

Species/Host

Rabbit

Reactivity

Equine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunohistochemical
staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human Thrombospondin 2. Please contact us for the exact immunogen sequence. The peptide is available as orb374869.

Target

Thrombospondin 2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

129 kDa

Uniprot ID
[Q03350](#), [P35440](#), [P35442](#)
NCBI
[21](#), [003247](#), [31](#), [003238](#)
Entrez
[7058](#)
Dilution Range

IHC-P: 1:200

Expiration Date

12 months from date of receipt.