

Product Datasheet

Integrin beta 1 antibody (orb251610)

Description


Integrin beta 1 antibody


Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Gallus, Human, Mouse, Porcine, Rat, Sheep


Conjugation

Unconjugated

Tested


ELISA, ICC, IF, IHC-P

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human Integrin beta 1. Please contact us for the exact immunogen sequence. The peptide is available as orb374870.


Target

Integrin beta 1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Immunohistochemical staining of pig stom...

Storage


Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG


Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Immunofluorescence image of pig stomach
...

MW

88 kDa

Uniprot ID
P05556, P09055, P49134
Entrez
3688
Dilution Range

IF/ICC: 1:200, IHC-P: 1:200

Expiration Date

12 months from date of receipt.