



Product Datasheet

VEGFR3 antibody (orb251613)


Description	VEGFR3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human VEGFR3. Please contact us for the exact immunogen sequence. The peptide is available as orb374872.
Target	VEGFR3
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	151 kDa
Uniprot ID	P35916, P35917, Q91ZT1
Entrez	2324
Dilution Range	IF/ICC: 1:200, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



IHC-P image of rat skin tissue using VEG...



IF analysis of mouse skin tissue using V...



IHC-P staining of rat skin tissue using ...