

Product Datasheet

BDNF antibody (orb251616)

Description

BDNF antibody




Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat



Conjugation

Unconjugated

Tested


ELISA, IHC-P

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human BDNF. Please contact us for the exact immunogen sequence. The peptide is available as orb374875.



Target

BDNF

Form/Appearance

10 mM PBS, 0.02% sodium azide


Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Immunohistochemical staining of paraffin...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.



Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P image of mouse skin tissue using a...

MW

28 kDa

Uniprot ID

[P23363](#), [P21237](#), [P23560](#)

NCBI

[001143805](#), [001137277](#), [11](#)

Entrez

[627](#)

Dilution Range

IHC-P: 1:200

Expiration Date

12 months from date of receipt.