



Product Datasheet

Iduronate 2 sulfatase antibody [Out of stock] (orb251621)


Description	Mouse monoclonal antibody to Iduronate 2 sulfatase
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-P, WB
Immunogen	Full length human recombinant protein of human IDS(NP_000193) produced in HEK293T cell.
Target	IDS
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol, 1% rAlbumin
Concentration	0.79 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Protein A/G purified
Entrez	3423
Dilution Range	FACS: 1:100, IF/ICC: 1:100, IHC-P: 10 µg/ml, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Br...



Immunofluorescence analysis of COS7 cell...



Western blot analysis of HEK293T cells u...