

## Product Datasheet

### Human EG-VEGF protein (orb251669)

|                        |                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Human EG-VEGF protein                                                                                          |
| <b>Reactivity</b>      | Human                                                                                                          |
| <b>Conjugation</b>     | Unconjugated                                                                                                   |
| <b>Target</b>          | Human EG-VEGF                                                                                                  |
| <b>Preservatives</b>   | Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 0.02 % Tween-20.                      |
| <b>Storage</b>         | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                                                          |
| <b>Purity</b>          | > 98 % by SDS-PAGE and HPLC analyses.                                                                          |
| <b>Source</b>          | E. coli.                                                                                                       |
| <b>MW</b>              | 9.7 kDa                                                                                                        |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                 |