

Product Datasheet

Rat MANF protein (orb251684)

Description	Recombinant Rat MANF protein
Reactivity	Rat
Conjugation	Unconjugated
Target	Rat MANF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 98 % by SDS-PAGE and HPLC analyses.
Source	E. coli.
MW	18.2 kDa
Expiration Date	6 months from date of receipt.