

Product Datasheet

Human UBE2M protein (orb251685)

Description	Recombinant Human UBE2M protein
Reactivity	Human
Conjugation	Unconjugated
Target	Human UBE2M
Preservatives	A 0.2 µm filtered concentrated solution in 50 mM HEPES, pH 6.5, with 125 mM NaCl, 10 % Glycerol, 5 % Trehalose, 1 mM DTT.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95 % by SDS-PAGE and HPLC analyses.
Source	E. coli.
MW	20.9 kDa
Expiration Date	6 months from date of receipt.