


Product Datasheet

NR2C2AP antibody (orb25172)


Description	NR2C2AP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Nuclear receptor 2C2-associated protein (7-139AA)
Target	NR2C2AP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86WQ0
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Immunohistochemical staining of human ov...



Western blot analysis of HL60 whole cell...