

Product Datasheet

PGLYRP4 antibody (FITC) (orb25324)

Description	PGLYRP4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Peptidoglycan recognition protein 4 protein (18-317AA)
Target	PGLYRP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96LB8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.