



Product Datasheet

ABCA2 Antibody (orb253650)

Description	ABCA2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BZC7
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of Human ce...



Immunohistochemical analysis of Human ga...