

Product Datasheet

POLR2K antibody (Biotin) (orb25386)

Description	POLR2K antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DNA-directed RNA polymerases I, II, and III subunit RPABC4 protein (1-58AA)
Target	POLR2K
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53803
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.