



Product Datasheet

PPP1CC antibody (orb25404)


Description	PPP1CC antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein phosphatase PP1-gamma catalytic subunit protein (151-323AA)
Target	PPP1CC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36873
Dilution Range	WB:1:2000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HEK293 whole c...



Immunohistochemistry of paraffin-embedde...