

Product Datasheet

PPP2CA antibody (Biotin) (orb25414)

Description	PPP2CA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform protein (1-309AA)
Target	PPP2CA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P67775
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.