

Product Datasheet

Prnp antibody (orb25434)

Description	Prnp antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rat Major prion protein (29-231AA)
Target	Prnp
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13852
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.