

Product Datasheet

Caspase 9 Mouse Monoclonal Antibody (orb254370)

Description	Caspase 9 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P55211
Expiration Date	12 months from date of receipt.