

Product Datasheet

alpha Lactalbumin Mouse Monoclonal Antibody (orb254377)

Description	alpha Lactalbumin Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P00709
Expiration Date	12 months from date of receipt.