

Product Datasheet

Oct1 Mouse Monoclonal Antibody (27-7G1) (orb254382)

Description	Oct1 Mouse Monoclonal Antibody (27-7G1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P14859
Expiration Date	12 months from date of receipt.