

Product Datasheet

Prss1 antibody (Biotin) (orb25448)

Description	Prss1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Anionic trypsin-1 protein (24-246AA)
Target	Prss1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00762
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.