



Product Datasheet

PRY antibody (orb25449)

Description	PRY antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human PTPN13-like protein, Y-linked protein (13-147AA)
Target	PRY
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14603
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue showing PRY antibody localization.



Immunohistochemical staining of human testis tissue showing PRY antibody localization.